


**BioPoly®**

# Radial Head



### Setting a New *Gold Standard* in Radial Head Replacement

The BioPoly Radial Head System offers 30 sizing options and features an industry-first Hydrophilic Polymer coupled with a smooth stem.

#### BioPoly Hydrophilic Polymer

##### Articulating Surface

- Ultra-High Molecular Weight Polyethylene (UHMWPE) and Hyaluronic Acid
- “Self-lubricates” in the joint due to its hydrophilic properties
- Interacts favorably with opposing cartilage while supporting physiological loads
- 85% lower Coefficient of Friction than Cobalt Chrome
- 40% less cartilage wear than Cobalt Chrome when articulating against cartilage

#### Polished Titanium Alloy

##### Smooth Stem


- Designed to self center during motion
- Freely aligns to patient anatomy
- Improved surgical efficiency






# BioPoly® Radial Head Product Catalog

## 30 Implant Sizing Options



10 Head + Stem combinations x 3 Height Offsets = 30 Available Sizes

## Reusable Instrument Set



## Implants

Description*	Part #
Radial Head Implant, 20 x 6(+0)mm	160-5060
Radial Head Implant, 20 x 6(+3)mm	160-5063
Radial Head Implant, 20 x 6(+6)mm	160-5066
Radial Head Implant, 20 x 7(+0)mm	160-5070
Radial Head Implant, 20 x 7(+3)mm	160-5073
Radial Head Implant, 20 x 7(+6)mm	160-5076
Radial Head Implant, 20 x 8(+0)mm	160-5080
Radial Head Implant, 20 x 8(+3)mm	160-5083
Radial Head Implant, 20 x 8(+6)mm	160-5086
Radial Head Implant, 23 x 6(+0)mm	160-5360
Radial Head Implant, 23 x 6(+3)mm	160-5363
Radial Head Implant, 23 x 6(+6)mm	160-5366
Radial Head Implant, 23 x 7(+0)mm	160-5370
Radial Head Implant, 23 x 7(+3)mm	160-5373
Radial Head Implant, 23 x 7(+6)mm	160-5376
Radial Head Implant, 23 x 8(+0)mm	160-5380
Radial Head Implant, 23 x 8(+3)mm	160-5383
Radial Head Implant, 23 x 8(+6)mm	160-5396
Radial Head Implant, 23 x 9(+0)mm	160-5390
Radial Head Implant, 23 x 9(+3)mm	160-5393
Radial Head Implant, 23 x 9(+6)mm	160-5396
Radial Head Implant, 26 x 7(+0)mm	160-5670
Radial Head Implant, 26 x 7(+3)mm	160-5673
Radial Head Implant, 26 x 7(+6)mm	160-5676
Radial Head Implant, 26 x 8(+0)mm	160-5680
Radial Head Implant, 26 x 8(+3)mm	160-5683
Radial Head Implant, 26 x 8(+6)mm	160-5686
Radial Head Implant, 26 x 9(+0)mm	160-5690
Radial Head Implant, 26 x 9(+3)mm	160-5693
Radial Head Implant, 26 x 9(+6)mm	160-5696

## Instruments

Description	Part #
Instrument Set	160-5500
Sizing Tray	160-2010-01
Instrument Handle	160-2010-02
Radius Retractor	160-2010-03
Broach - 6mm	160-2011-06
Broach - 7mm	160-2011-07
Broach - 8mm	160-2011-08
Broach - 9mm	160-2011-09
Planer - 6/7mm	160-2012-06
Planer - 8/9mm	160-2012-08
Height Trial	160-2013-01
Stem Trial - 6mm	160-2014-06
Stem Trial - 7mm	160-2014-07
Stem Trial - 8mm	160-2014-08
Stem Trial - 9mm	160-2014-09
Head Trial - 20mm	160-2015-20
Head Trial - 23mm	160-2015-23
Head Trial - 26mm	160-2015-26
Offset Trial - +3mm	160-2016-03
Offset Trial - +6mm	160-2016-06
Tamp	160-2017-01

\*[Head Ø x Stem Ø (offset)]



**Advancing Materials. Advancing Outcomes.™**

Manufactured By:

BioPoly LLC, a Schwartz Biomedical Company  
7136 Gettysburg Pike, Fort Wayne, IN 46804, USA  
Phone 260.999.6135 | [www.BioPolyortho.com](http://www.BioPolyortho.com)

This product is covered by US Patent No. 7,662,954; 10,405,982; 11,690,722 and other patents pending. Copyright 2024 BioPoly LLC. All rights reserved. BioPoly is a trademark of BioPoly LLC. Printed in the USA.